




PRECISION ON GLASS

IMT Glass Microfluidics for Life Science and Diagnostics


Flexible process offerings that enable microfluidic solutions in glass: design consultancy, prototyping, and scalable manufacturing

More informations?
Contact us:
info@imtag.ch







Nanopatterning enabling high-throughput multiplex assays and next generation sequencing



Structured waveguides enabling signal-to-noise reduction



Complex glass structuring enabling multiple-emulsion droplet generation



Electrode integration enabling pathogen detection

Applications

NGS flow cell, organ-on-a-chip, lab-on-a-chip, single-cell analysis, cell enrichment, sample preparation and many more

Advantages of microfluidics in glass

- Selection of available glass types
- Bioinert
- Excellent chemical, mechanical, and optical properties
- Outstanding surface properties
- Cost-efficient scaling from prototyping to mass manufacturing

Processes

- Nano-scale patterning of glass down to 350 nm
 - Pillars, wells, channels
- Electrode integration
 - Materials: Au, Pt, ITO, Ti
 - Features sizes down to 2 μm
- Biocompatible bonding
- Structured Bio-functionalisation
- ISO 9001 : 2015

About IMT

- Foundry for optical, electrical and microfluidic structures and components
- Fast prototyping through in-house mask manufacture
- Staff of 110 employees
- 1600 m² clean room

IMT

IMT Masken und Teilungen AG
 Im Langacher 46 · 8606 Greifensee (Switzerland)
 T +41 44 943 19 00 · info@imtag.ch · imtag.ch